

Eleganza Classic



Hospital Bedside Cabinet

Classic - universal equipment for hospital rooms

The Eleganza Classic patient cabinet with or without collapsible dining board allows placement on the left or right side of the bed while maintaining all its functions.

- HEIGHT ADJUSTMENT AND TILTING OF THE EATING BOARD
- EASY MANIPULATION DUE TO THE DOUBLE CASTORS
- EASY CLEANING



DOUBLE BEDSIDE USE

- The cupboards are accessible from the front and back of the cabinet and drawers can slide out in both directions
- The mechanism of the dining board is integrated into the body
- Upper drawer, niche, lower space with sliding shelf and doors.
- The cabinets are equipped with 4 double castors (2 include a brake)
- Models without eating board also feature integrated holder for plastic bottles and magazines.

ADJUSTMENT

- The mechanical system of the eating board is directly integrated in the corpus of the cabinet.
- Painted galvanized metal cabinet body which guarantees high resistance to wear, corrosion and peeling. The top and eating boards are made from high pressure laminate, with high mechanical, chemical and thermal durability.
- The top boards of the cabinet feature a gallery, which also serves as a handrail for the positioning of the table. Both boards have raised edges that prevent object from falling down.



Integrated plastic holder for three bottles for model without eating board.



The plastic water resistant organisers can be taken out from metal holding frame.



Top board with raised edges and handrails for manipulation. Side towel holder.

TECHNICAL PARAMETERS

Outer dimensions	
(without eating board)	50 x 51,5 cm
Eating board dimensions	
(between edges)	60 x 34,5 cm
Table height (upper desk)	86,5 cm
Height of eating board	
(castors diam. 7,5 cm)	72 – 108 cm
Max. load of eating board	
(in the centre)	9 kg
Castors	7,5 cm
Weight	35 kg