

## Product specification:

### Patient's bedside tables with wooden finishing

model No: 178A8118110A1 and 47528118D0L01



#### Description

#### Technical parameters

#### model No: 178A8118110A1

- drawers and cabinets accessible from both sides of locker
- eating board (height adjustment, tilt)
- made from melamin plated wooden boards
- standard colour – beech 01
- metal parts RAL 9006 – silver grey
- in quantity of 350 pcs can be offered in different colour – please specify (send a sample)
- also can be equipped by locks – specify if you want to lock shelf (1 lock) or the cabinet (2 locks needed)

total dimensions (w x d x h)	51 x 48 x 89 cm
eating board dimensions	60 x 36 cm
height of eating board	82 – 110 cm

#### model No: 47528118D0L01

- drawers and cabinets accessible from one side (left/right - 47528118D0P01) of locker
- made from melamin plated wooden boards
- standard colour – beech 01
- in quantity of 350 pcs can be offered in different colour – please specify (send a sample)
- also can be equipped by locks – specify if you want to lock shelf (1 lock) or the cabinet (2 locks needed)

total dimensions (w x d x h)	47x 40x 83cm
------------------------------	--------------